

Some notes about scheduler functions:

The top of list shows what is currently executing - it is a FIFO stack except for the default experiment which always sits at the top and pauses for other experiments that execute. While experiments are waiting to execute in the stack (0%). To determine if an experiment is pauseable, a yeild command is specifically included in the experiment script.

There are other experiments that are not in the stack -- they execute every XX seconds. When they count down, they break into the current experiment if it is pauseable, else they jump to the top of the stack.

Scheduler	
Name	Status
Default	PAUSED
PhaseFlopBatch	10%
PhaseFlopRamsey	82%
CoolingExpt	0%
SqueezingExpt	0%
Field Set	96s
Mode Check	539s

**Buttons:**

- Remove
- Pause
- Restart
- Find
- Abort

Find Button:  
Selects and  
displays Expt.  
Explorer window  
for a given  
scheduled  
experiment

Abort Button: clears  
the queue and  
schedules the  
experiment stored  
in the "abort"  
parameter in the  
database.

## Experiment List

Debugging  
Rand. Bench  
ODF  
Transport  
+

- + 1 Be
  - Expt One
  - Expt Two
  - Expt Three
- + 2 Mg
  - MS Gate
  - Sz Gate
  - CZ Gate
  - I Gate
- + Hardware
  - TTL
  - DDS
  - Waveplate
  - etc.

Run

Schedule

Parse

Automatically Parse

A menu that opens when RightClicking on an experiment:  
Show: open Expt Explorer without parsing.  
Duplicate: make a copy of experiment  
Move: Move experiment to different tab  
Delete: delete experiment

Selecting an experiment will generate the GUI associates with the experiment.

These experiment GUI (widgets) are highly configurable for each experiment: which widgets to be prompt and what are the functionality of those widgets.

These widget can be:

Sweep control  
Plot control  
Fit control

Run: Will immediately start running selected experiment if currently running is pauseable, else will go to top of the queue.

Schedule: Add to the bottom of the queue -- but will run before experiments that are waiting on a cycle time.

Parse: Will parse the experiment and show the Experiment Explorer window

Parameter database contains all the information needed to run experiments.

Entries can be defined as more than just scalars, such that the database can store reference histograms, PDQ waveforms, density matrices, matrices of parameters (e.g. the spin flip frequency each ion in a multi-ion experiment, or a matrix of the coupling parameters in a multi-ion experiment).

Parameter Database		
<input type="checkbox"/> General		
B field	0.1234	
TempTable1	20.1	
TempTable2	19.8	
<input checked="" type="checkbox"/> Rotating Wall		
<input type="checkbox"/> Transport		
Wvf 1.1	dtype[vect]	
Wvf 3.4	dtype[vect]	
<input type="checkbox"/> Reference Histograms		
BeHist	dtype[vect]	
<input type="checkbox"/> Tomography		
<input type="checkbox"/> 2 Mg SigmaZ Gate		
Density Matrix	dtype[matr]	
Histogram	dtype[vect]	

search

Config

Config: Allow changing of the properties of the selected parameter:  
name  
class  
numeric format: 0.01 or 1x10^-2  
Console update configuration: on/off, frequency

Prints to the console have priority, and the console can be figured to show windows of different level of priority. So in one window, only high priority messages that you want to be able to be easy to read will show -- errors and fits for examples.

Can choose color of prints to the console

Other windows show full console information

The screenshot shows a window titled "Pytheas Console" with a black background and white text. It displays several types of messages:

- ERROR!** Experiment parse failed because of error code 9909!
- ERROR!** Other error messages I don't know about.
- FIT RESULTS:** Fit Routine Quadratic:  $y(x) = K_0 + K_1*x + K_2*x^2$
- $K_0 = 0 \pm \dots$
- $K_1 = 1.2 \pm 0.1$
- $K_2 = 2.0 \pm 0.01$
- 10274
- 1209
- 14.65670
- Print something
- update something
- other print
- print print print
- ERROR!** Experiment parse failed because of error code 9909!
- other print
- 14.65670
- Print something
- ERROR!** Other error messages I don't know about.
- Print something
- update something
- other print
- FIT RESULTS:** Fit Routine Quadratic:  $y(x) = K_0 + K_1*x + K_2*x^2$
- $K_0 = 0 \pm \dots$
- $K_1 = 1.2 \pm 0.1$
- $K_2 = 2.0 \pm 0.01$

At the bottom of the window is a search bar with the placeholder text "search".

How hard is search to implement?  
Useful? Low priority

[Experiment](#)

[Parameter](#)

[Database](#)

[Console](#)

[About](#)

X

Scheduler

[Experiment List](#)

Prompt a window asking:

"Are you sure you want to quit?"

and gives you the option to save the layout.

\*\* Other window X buttons will NOT close  
the GUI, just the window \*\*


Your Fav Code Editor

```
class scan1Dexp():
    parameters = "scanned_experiment scanned_parameter"

    def get_widgets():
        return [
            Checkbox("2nd Doppler Cooling", "enable_doppler_cool2"),
            GlorifiedSlider(self.scanned_parameter) ] +
    self.scanned_experiment.get_widgets()

    def run(self):
        for x in range(start, stop)
            exp = self.scanned_experiment(self, **{self.scanned_parameter: x})
            yield exp
            broadcast(s, self.histogram)
            plot(xdata, ydata, o)
        if enableSave:
            nameAppend = getTime()
            self.saveData(ParamDB.saveFile+nameAppend)
```


SB has said that a class like scan1Dexp() would likely be something in a library provided by ARTIQ, so many of the details of how this would work would not need to be accessed day to day. Still, it shows the simplicity of adding widgets as your own classes can have a get\_widgets() function.



## General Description


The master control panel of the GUI is a tab bar/menu bar with a list of items in different categories.

A separate window will open when click on these items.



## Pytheas GUI Scheduler

Some notes about scheduler functions:  
The top of list shows what is currently executing - it is a FIFO stack except for the default experiment which always sits at the top and pauses for other experiments that execute. While experiments are waiting to execute in the stack (0%), they can be reorganized (^ v buttons beside them). There are other experiments that are not in the stack -- they execute every XX seconds. When they count down, they break into the current experiment if it is pauseable, else they jump to the top of the stack.



A menu that opens when RightClicking on an experiment:  
Show: open Expt Explorer without parsing.  
Duplicate: make a copy of experiment  
Move: Move experiment to different tab  
Delete: delete experiment

The experiment control panel will change to the corresponding selected experiment in this list.

Ting Rei's note: we should also include the option that a separate experiment control panel will be shown instead of the "default experiment control panel"

Ting Rei's comments:  
1. how to implement automatic loading?  
2. how to implement "intertwine" different experiment

JGB: I'm not sure what you mean by automatic loading.

Ting Rei: justin, this is a note I wrote to remind myself. You can ignore it if you want.

This first window is the main way of viewing and interacting with the parameter database. Entries can be defined as more than just scalars, such that the database can store reference histograms, PDQ waveforms, density matrices, matrices of parameters (e.g. the spin flip frequency each ion in a multi-ion experiment, or a matrix of the coupling parameters in a multi-ion experiment).

Reference histogram is saved every xxx minutes. User can access and look at how the each reference histogram behave.

History Plot: A simple plot of the selected parameter vs time. Can be used as a "logging" type plot. Of course, this database is on Git, so I suppose it can be accessed and plotted through other means (grafana? others?) but this could provide a simple interface built in. -- How to handle other datatypes?

History List: Display a .CSV of last 10 (20? 50?) data points for a parameter

Export: A quick way to extract data from the database and save it to another file

**Pytheas GUI Parameter Database**

Category	Parameter	Type
General	B field	0.1234
	TempTable1	20.1
	TempTable2	19.8
Rotating Wall	Wvf 1.1	dtype[vect]
	Wvf 3.4	dtype[vect]
Transport	BeHist	dtype[vect]
	Density Matrix	dtype[matr]
Tomography	Histogram	dtype[vect]
	2 Mg SigmaZ Gate	

Buttons: History Plot, History List, Export..., Config

Search: search

Config: Allow changing of the properties of the selected parameter:  
name  
class  
numeric format: 0.01 or 1x10^-2  
live feed configuration: on/off, frequency

These are 2 other windows connected to the database, but the idea is that they are "Display Only" to reduce complexity associated with interacting with the database.

A live feed of parameters that have been updated, regardless of which experiment or which protocol that update them.

I removed the buttons here with the idea that this is just a 'display' -- to interact with a parameter you have to go to the database window.

**Pytheas GUI Database History**

Parameter	Value	Last Update
2110 CO AO Freq	265	14:54:33:265
<b>B Field</b>	<b>1.265</b>	<b>14:53:23:100</b>
Trap RF amplitude	25	14:49:00:298
Doppler Waveplate	245.1	14:46:10:987
etc.		
etc.		

showing the last 30 (or any number) of parameters that was last changed

Description: Show the description of the selected parameter above.  
E.g. which experiment change this parameter and what fitting routine was used.

**Pytheas GUI Expt Parameters**

myExperimentName

Parameter	Value
2110 Co AO Freq	265
2110 90 AO Phase	192
Separate Waveform 1	wvf 1.1
<b>Doppler Waveplate</b>	<b>245.1</b>

Pytheas GUI Experiment Explorer

# myExperimentName

Pauseable? Repeat every 60 sec Max Pause 40 Num Expts 40 Build

Experiment definition file  
 Run C:/qc/exptFiles/SidebandCooling.  
 Save? Output: C:/qc/userData/output/

Expt Variables ▲ Value Scan-able?  
 arm\_t 10.00 Y  
 final\_p 0 Y  
 pi\_time 32.447066 Y  
 AO pulse shaping H N  
 Be Separate Wvf Wvf 1.1 N

Profile Select: DetectionShelving  
 Sideband Cooling  
 Doppler Only

Enable Scan   
 Variable to scan: arm\_t

Start Stop Points Scans Random?  
 0 15000 6 1

<-> CenterOn Cursor ->-

Scan history  
 Parameter▲ Value◀ Last Scan  
 tFlop 15.6 Today 09:35:10  
 Freq 254 Today 09:40:20  
 LabBrick 1 -0.265 Yesterday 12:20:20

Data Linker Plotting Debug

QuickPlot

X Data	Y Data	Histograms
Scan Variable	Scan Variable	PMT1 []
PMT1_avg	PMT1_avg	PMT2 []
PMT1_std	PMT1_std	
NI_AIO2_time	NI_AIO2_time	
NI_AIO2_avg	NI_AIO2_avg	

Fit function: Choose Y data to fit: ComboBox  
 Choose fit function:  $y(x) = x^2$   
 Refit  
 $y(x) = \sin(x)$

Param	Guess	Fitted	Fitted Err	Fixed?
K0	0	0	0	<input checked="" type="radio"/>
K1	1	1.2	0.05	<input type="radio"/>

Analysis file   
 Filename: C:/qc/exptFiles/SidebandCooling.dc Run

Data Linker Plotting Debug

A tab for anything related to debugging the sequence of code for the experiment. More ideas for useful GUI things?

Complie Error list  
 Item One  
 Item Two  
 Item Three

Num Errors 3  
 Num Pulses 40

Pulse Visualizer  
 Visualize  Visualize Depth 1  
 Just a place holder graphic

Plotting Tab:

Idea is that Quick plot is only for making 1 simple X vs Y plot (and displaying histograms). Make this DEAD simple. If more complex analysis is to be done, open up the text file below.

In the code, have to define certain variables as data that can be plotted, then it shows up on list.

Build is used if edits to the file would cause a change to the GUI, for example when the experiment is first created or if some new kind of tab is added.

Visualize pops open a window to show the experiment pulses. Depth specifies complexity of the graph. Could use flags built into pulses in the experiment code to provide the different levels.


# TTL Switch Panel



## TTL Grouping

Be Raman All  
Be Doppler All  
Mg Co  
ODF  
etc.

Add    Configure    Clear



Clicking the button to rotationally force on, off or expt. control. The color on the button show the status of that TTL

The blue/purple indicators shows the current status of the TTL. They are update in real time.

All On    All Off    Expt. Control

Configure

Description of the TTL (User define)

Change description, how the status indicator shows the status of the TTL, etc...

## DDS Control Panel

### DDS 0 (Be Doppler)

Frequency 250

Amplitude -0.5 dBm

Phase 90

Reset

Manual Override

### DDS 1 (Mg Co)

Frequency 327

Amplitude 3.2 dBm

Phase 0

Reset

Manual Override

Showing the DDS's value that is being set to

Update the value at the end of each experiment: the idea is that we can see this value is changing in an experiment where a frequency/phase is being scanned.

User can change the DDS values when the "Manual Override" check box is checked.

When checked, user can manually override the DDS setting even when the experiment that is running try to change the setting

When left uncheck, the properties is controlled by the running experiment

### Master Control

Frequency 200

Amplitude 0 dBm


Phase 0

Reset All

Manual Override

Configure

} Control the setting of all DDS



Clicking on the configure button will pop up the experiment control windows which is similar to the "default" experiment control window, except there will be an obvious indication that this experiment window is setting the experiment control of the Bootup Experiment, the Idle experiments or the Power failure experiment.

In principle, one can even group every "special condition" experiment into a selection of experiment with different given condition to run. The windows on the left will enforce that only one experiment will run in the situation of bootup and power failure.

Pytheas Hardware Dashboard

**TTL myname**

**TTL Grouping**

- Be Co (TTL 0)
- Be 90 (TTL 1)
- uWv 1 (TTL 2)
- BD (TTL 3)
- Mg 90 (TTL 4)
- ODF (TTL 5)
- TTL 6
- TTL 7

**Add** **Configure**

**All On** **All Off** **Configure**

**NI\_AIO1**

**Hardware Config File** /NI\_artiqDriver

**Clock Select** External (Pin 1)

**AO Outputs**

Enable?

Chn	ProgName	Value
0	DC1	0.000
1	DC2	0.000
2	DC3	1.000

**AO Inputs**

Enable?

Chn	Dataname
0	HetPhase1
1	Bfield1

Instead of an "output" list and an input list, I think everyone piece of hardware (or network driver) should get 1 entry. The indicator tells if it is connected without an error. Then in its dashboard you define whether it is an input, an output, or both!

Buttons to shelf back into the menu, or pop out to a separate window.

The "input" is list of data collected by the ARTIQ such that we can look at their behavior when we scan certain variable in experiment.

The general idea is: ARTIQ does not only output PMT counts, but one can use ARTIQ to look at other parameters. Moreover, one can record the behavior of 2 parameters as an user scans a variable in an experiment, this allows us to look at correlation.

Buttons to shelf back into the menu, or pop out to a separate window.

**Pytheas GUI Experiment Control**

### myExperimentName

Experiment definition file

**Scan Variable**

Expt Variables	Value
arm_t	10.00
final_p	32.447066
initial_p	
middle_p	
pi_time	

Many scan variable in this list will be overlapping with the expt variables.

This box allow user to set choose which variable is scanable and which can only be changed but not scanable.

**Enable Scan**

Number of Scan: 1

Start Stop Interval Random?

0 15000 0.01

<-> CenterOn Cursor ->->

**Scan history**

Parameter	Value	Last Scan
tFlop	15.6	Today 09:35:10
Freq	254	Today 09:40:20
LabBrick 1	-0.265	Yesterday 12:20:20

Simple 2D Plot Roll Plot

If Roll Plot Roll Plot Max time

As a fail safe if something will go wrong at long time.

**Measurement X Plotting X Quick Control X Repeat X Saving X Analysis X**

**QuickPlot**

Fit Type: Sine Curve (simple)

Lorentzian Exponential Parabola

a Sin(x+phi) +c

Show the functional form of the fit function

Parameter Value Error

a	0.86	0.002
phi	0.25	0.006
c	0.652	0.02

Ting Rei's notes:  
Maybe another window for plotting control and analysis?

The "Analysis" tab is for using customized analysis routine (Python or other software) to analyze the data.

**Measurement Plotting Quick Control Repeat Saving Analysis Alerting**

**PMT 1**

- Fluorescence C:/qc/Out/SidebandCooling/PMT1\_Fluorescence.csv Setting
- Histogram C:/qc/Out/SidebandCooling/PMT1\_Histogram.csv Setting
- 2 Be Reference Histogram Maximum Likelihood Fit C:/qc/Out/SidebandCooling/PMT1\_RefHistFit.csv Setting
- Agnostic Analysis 1 C:/qc/Out/SidebandCooling/PMT1\_Agnostic1.csv Setting

**ADC Channel 1**

- Temperature C:/qc/Out/SidebandCooling/ADC1\_temp.csv Setting

**Measurement Plotting Quick Control Repeat Saving Analysis Alerting**

PMT 1
 

- Show Counts
- Show Histogram
- Update to reference histogram

PMT 2
 

- Show Counts
- Show Histogram
- Update to reference histogram

ADC Channel 1 (AOM temperature)
 

- Show Value

ADC Channel 2 (Room temperature )
 

- Setting

This is a list of measurement input that can be selected as the plotting (and data saving) dependent variables.

If more than one variable is selected, two plotting graph will pop up as a variable is being scanned.

**Measurement Plotting Quick Control Repeat Saving Analysis Alerting**

Pausable

Repeat 300s Setting

Automatic Loading if loss ion Setting

Rescue Ion Setting

This panel shouldn't be called "Repeat", but I don't really have a good name for it right now.

The function of this panel is to put extra feature to the experiment, and these features can be add/remove by user depending on which experiment.

For example: these features can be:

1. "Rescue ion": more for detuned Doppler cooling if the fluorescence is below certain threshold
2. "Automatic loading": Load ion automatically in the event of loss ion (this will be very useful for long algorithm experiment such as RB)
- 3.

**Measurement Plotting Quick Control Repeat Saving Analysis Alerting**

**Alert me if**

- PMT 1 Fluorescence Count below 20
- Fit Fail
- ADC Channel 7 below 0.85
- TTL Input 4 is High etc.

Alert Method:

- Sticky Note on Console
- Console Update
- Prompted Window
- TTL 9 output HIGH
- Record in Data Save File
- Email Me tingrei.tan@nist.gov
- Text Me
- Clash the GUI
- etc.

**Configuration**

# myExperimentName

Experiment definition file

**Run**

C:/qc/exptFiles/SidebandCooling

 Save? Output: C:/qc/userData/output

Expt Variables	Value
arm_t	10.00
pi_time	32.447066
Freq	236.23
Shim 38	-0.0365
AO Pulse Shaping	Enable
Be Separate Wvf	Wvf v1.1

## Scan Variable 1

Start	Stop	freq	Random?
0	15000	pi_time	<input type="checkbox"/>
<input type="button" value="&lt;--&gt;"/> <span style="margin: 0 10px;">CenterOn Cursor</span> <input type="button" value="-&gt;&lt;-"/>			

## Scan Variable 2

Start	Stop	Interval	Random?
0	15000	0.01	<input type="checkbox"/>
<input type="button" value="&lt;--&gt;"/> <span style="margin: 0 10px;">CenterOn Cursor</span> <input type="button" value="-&gt;&lt;-"/>			

## Scan history

Parameter	Value	Last Scan
tFlop	15.6	Today 09:35:10
Freq	254	Today 09:40:20
LabBrick 1	-0.265	Yesterday 12:20:20

Simple

If Roll Plot

2D Plot

Roll Plot Max time

Roll Plot

Number of Scan:

1

Measurement X

Plotting X

Quick Control X

Repeat X

Saving X

Analysis X

## QuickPlot

Fit Type

Sine Curve (simple)

Lorentzian

Exponential

Parabola

 $a \sin(x+\phi) + c$ 

Show the functional form of the fit function

Parameter

Parameter	Value	Error
a	0.86	0.002
phi	0.25	0.006
c	0.652	0.02

# myExperimentName

Experiment definition file

C:/qc/exptFiles/SidebandCooling.dc

Save? Output: C:/qc/userData/output/